



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD107b Monoclonal Antibody

单克隆抗体

产品货号: MA02244

第三版

免疫原:Fusion protein of LAMP2

亚型:IgA

克隆:34

描述:Lysosomal-associated membrane protein 2 (LAMP2,synonyms: LAMPB,CD107b) is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. LAMP2 may plays a role in tumor cell metastasis. It may also functions in the protection,maintenance,and adhesion of the lysosome. Prior to posttranslational modification,LAMP2 is a ~45 kDa polypeptide. Mature,functional LAMP2 is extensively glycosylated with a variety of different N linked and O linked oligosaccharides.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF