

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MAOB Monoclonal Antibody

单克隆抗体

产品货号: MA02247

第三版

免疫原:Fusion protein of MAOB

亚型:IgG1

克隆:31

描述:The MAOB gene encodes a 520-amino acid protein with a molecular mass of 58.8 kD and shows 70% amino acid identity to MAOA. MAOA and MAOB are present in the outer mitochondrial membrane in the central nervous system and peripheral tissues. A comparison of highly purified human placental MAOA and human liver MAOB revealed that the A form of the enzyme is larger by 2 kDa and only one potential N-glycosylation site exists in each protein, with each site in a different relevant position.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF