



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse MAPK13 Monoclonal Antibody

单克隆抗体

产品货号: MA02248

第三版

免疫原:Fusion protein of MAPK13

亚型:IgG1

克隆:224

描述:MAPK13(Mitogen-activated protein kinase 13) is also named as PRKM13,SAPK4 and belongs to the protein kinase superfamily. MAPK13 is one of the four p38 MAPKs which play an important role in the cascades of cellular responses evoked by extracellular stimuli such as proinflammatory cytokines or physical stress leading to direct activation of transcription factors such as ELK1 and ATF2.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF