

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MBIP Monoclonal Antibody

单克隆抗体

产品货号: MA02249

第三版

免疫原:Fusion protein of MBIP

亚型:IgG2b

克隆:126

描述:MBIP (MAP3K12 binding inhibitory protein 1),also known as MAPK upstream kinase (MUK)-binding inhibitory protein,is a 344 amino acid nuclear and cytoplasmic protein that assists in the activation of the JNK pathway by inhibiting ZPK activity. Ubiquitously expressed,MBIP is found at high levels in lung and heart and is one component of the ADA2A-containing complex (ATAC) in conjunction with ADA3,GCN5,CDC101,TADA2L,WDR5,YEATS2,DR1 and CRP2BP.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF