

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MYOD1 Monoclonal Antibody

单克隆抗体

产品货号: MA02253

第三版

免疫原:Fusion protein of MYOD1

亚型:IgG1

克隆:654

描述:This gene encodes a nuclear protein that belongs to the basic helix-loop-helix family of transcription factors and the myogenic factors subfamily. It regulates muscle cell differentiation by inducing cell cycle arrest, a prerequisite for myogenic initiation. The protein is also involved in muscle regeneration. It activates its own transcription which may stabilize commitment to myogenesis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF