

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NENF Monoclonal Antibody

单克隆抗体

产品货号: MA02254

第三版

免疫原:Fusion protein of NENF

亚型:IgG2a

克隆:6540

描述:This gene encodes a neurotrophic factor that may play a role in neuron differentiation and development. A pseudogene of this gene is found on chromosome 12. Alternate splicing of this gene results in multiple transcript variants.NENF (Neudesin Neurotrophic Factor) is a Protein Coding gene. GO annotations related to this gene include growth factor activity and heme binding. An important paralog of this gene is CYB5D2.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF