

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GNAS Monoclonal Antibody

单克隆抗体

产品货号: MA02255

第三版

免疫原:Fusion protein of NESP55/GNAS

亚型:IgG2b

克隆:364

**描述:**Guanine nucleotide binding protein (G protein),alpha stimulating activity polypeptide 1 (GNAS1) is the ubiquitously expressed heterotrimeric G protein that couples receptors to the effector enzyme adenylyl cyclase and is required for receptor-stimulated intracellular cAMP generation. Mutations of Gs(alpha) residues involved in the GTPase reaction that lead to constitutive activation are present in endocrine tumors, fibrous dysplasia of bone, and McCune-Albright syndrome. The molecular weight of Gs(alpha) protein is about 46 kDa.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF