

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat OGDH Monoclonal Antibody

单克隆抗体

产品货号: MA02257

第三版

免疫原:Fusion protein of OGDH

亚型:IgG2b

克隆:160

描述:OGDH(2-oxoglutarate dehydrogenase,mitochondrial) is also named as OGDC-E1and belongs to the alpha-ketoglutarate dehydrogenase family. It is one of the three enzymatic components of 2-oxoglutarate dehydrogenase complex,which catalyzes the overall conversion of 2-oxoglutarate to succinyl-CoA and CO₂. This protein has 3 isoforms produced by alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF