



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human OXA1L Monoclonal Antibody

单克隆抗体

产品货号: MA02258

第三版

免疫原:Fusion protein of OXA1L

亚型:IgG1

克隆:320

描述:OXA1L/Oxa1 is a member of the conserved Oxa1/YidC/Alb3 protein family and is required for cytochrome c oxidase assembly. The OXA1L protein was initially identified in yeast proteins involved in the biogenesis of the mitochondrial respiratory chain complexes. It locates in the inner membrane of mitochondria and may be involved in mediating the insertion of proteins from the matrix into the inner membrane.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF