



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PA2G4 Monoclonal Antibody

单克隆抗体

产品货号: MA02259

第三版

免疫原:Fusion protein of PA2G4

亚型:IgG2a

克隆:9544

描述:PA2G4 (EBP1),first identified as an ErbB3 binding protein,is a 38 kDa protein that is widely expressed in cultured cells and tissues. The presence of several functional domains such as a nuclear localization sequence (NLS),LxxLL and LxCxE suggests the involvement of Ebp1 in cell signaling pathways and gene transcription regulation. This antibody is a mouse mono-clonal antibody raised against full length PA2G4 of human origin.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF