

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SERPINE1 Monoclonal Antibody

单克隆抗体

产品货号: MA02260

第三版

免疫原:Fusion protein of PAI-1

亚型:IgG1

克隆:

描述:SERPINE1,also named as Plasminogen activator inhibitor type 1 (PAI-1),is a member of the serine protease inhibitor (SERPIN) superfamily. PAI-1 is the principal inhibitor of tissue plasminogen activator (tPA) and urokinase (uPA),and hence is an inhibitor of fibrinolysis. Defects in this gene are the cause of plasminogen activator inhibitor-1 deficiency (PAI-1 deficiency),and high concentrations of the gene product are associated with thrombophilia.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF