

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SDCCAG8 Monoclonal Antibody

单克隆抗体

产品货号: MA02269

第三版

免疫原:Fusion protein of SDCCAG8

亚型:IgG2b

克隆:447

描述:SDCCAG8,also named as CCCAP and NY-CO-8,plays a role in the establishment of cell polarity and epithelial lumen formation. It may play a role in ciliogenesis. Loss of SDCCAG8 function as a cause of a retinal-renal ciliopathy and validates exome capture analysis for broadly heterogeneous single-gene disorders. SDCCAG8 is localized at both centrioles and interacts directly with OFD1 (oral-facial-digital syndrome 1),which is associated with NPHP-RC.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF