



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat SERPINA10 Monoclonal Antibody

单克隆抗体

产品货号: MA02270

第三版

免疫原:Fusion protein of SERPINA10

亚型:IgG2a

克隆:156

描述:SERPINA10 belongs to the serpin family. It is predominantly expressed in the liver and secreted in plasma. It inhibits the activity of coagulation factors Xa and XIa in the presence of protein Z, calcium and phospholipid. Mutations in this gene are associated with venous thrombosis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF