

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SGTA Monoclonal Antibody

单克隆抗体

产品货号: MA02272

第三版

免疫原:Fusion protein of SGTA

亚型:IgG2b

克隆:

描述:Small glutamine-rich tetratricopeptide repeat-containing protein alpha (SGTA) is a steroid receptor molecular co-chaperone critical in a range of biological processes,including viral infection,cell division,mitosis,and cell cycle checkpoint activation. Overexpression of SGTA has been found in various cancers,including prostate/ovary/testis/liver/lung cancer.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF