

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ST6GALNAC6 Monoclonal Antibody

单克隆抗体

产品货号: MA02273

第三版

免疫原:Fusion protein of ST6GALNAC6

亚型:IgG1

克隆:186

描述:ST6GALNAC6 is a sialyltransferase which is responsible for the synthesis of all a-series gangliosides. It has a broad expression pattern. ST6GALNAC6 was reported to be a key enzyme in the synthesis of disialyl monosialyl Lewis (Lea), a representative tumor-associated antigen in cancers of the pancreas and colon.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF