

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SULT1A1 Monoclonal Antibody

单克隆抗体

产品货号: MA02275

第三版

免疫原:Fusion protein of SULT1A1

亚型:IgG2a

克隆:28

描述:SULT1A1,also named as ST1A1,P-PST 1,ST1A3,Ts-PST,P-PST 1,STP and STP1,belongs to the sulfotransferase 1 family. It catalyzes the sulfate conjugation of catecholamines,phenolic drugs and neurotransmitters. SULT1A1 is also responsible for the sulfation and activation of minoxidil. It mediates the metabolic activation of carcinogenic N-hydroxyarylamines to DNA binding products and could so participate as modulating factor of cancer risk. SULT1A1 is one of two phenol sulfotransferases with thermostable enzyme activity.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF