

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SULT4A1 Monoclonal Antibody

单克隆抗体

产品货号: MA02276

第三版

免疫原:Fusion protein of SULT4A1

亚型:IgG1

克隆:

描述:SULT4A1,also named as SULTX3,NST,ST4A1,hBR-STL and hBR-STL-1,belongs to the sulfotransferase 1 family. It may catalyze the sulfate conjugation of many drugs,xenobiotic compounds,hormones,and neurotransmitters. SULT4A1 displays activity towards L-triiodothyronine,thyroxine,estrone,p-nitrophenol,2-naphtylamine,and 2-naphtol. It may have a role in the sulfation of drugs and neurotransmitters in the CNS. SULT4A1 is a brain-specific sulfotransferase believed to be involved in the metabolism of neurotransmitters. It is a novel cytosolic sulfotransferase that is primarily expressed in the brain.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF