

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TRAPPC9 Monoclonal Antibody

单克隆抗体

产品货号: MA02279

第三版

免疫原:Fusion protein of TRAPPC9/NIBP

亚型:IgG2b

克隆:538

描述:TRAPPC9,also named as KIAA1882 and NIBP,belongs to the NIBP family. It functions as an activator of NF-kappa-B through increased phosphorylation of the IKK complex. TRAPPC9 may function in neuronal cells differentiation and play a role in vesicular transport from endoplasmic reticulum to Golgi. TRAPPC9 was found in neurons of the cerebral cortex,hippocampus,and deep gray matter.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF