



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TRIM44 Monoclonal Antibody

单克隆抗体

产品货号: MA02280

第三版

免疫原:Fusion protein of TRIM44

亚型:IgG2b

克隆:573

描述:TRIM44 is one of the family member of TRIM protein,which contains an N-terminal ubiquitin hydrolase-type zinc-finger domain, followed by a coiled-coil domain, a zinc-finger B-box homology domain, and a second coiled-coil domain near the C-terminus. Members of the TRIM protein family are involved in various cellular processes, such as cell proliferation, oncogenesis and antiviral defense.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF