



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TRIP6 Monoclonal Antibody

单克隆抗体

产品货号: MA02281

第三版

免疫原:Fusion protein of TRIP6

亚型:IgG1

克隆:598

描述:TRIP6 (thyroid receptor-interacting protein 6),also known as ZRP-1 (zyxin-related protein 1),is a member of the zyxin family that has been implicated in actin reorganization and cell motility. Trip6 has been proposed to transport signals from the cell surface to the nucleus.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF