

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human USP5 Monoclonal Antibody

单克隆抗体

产品货号: MA02283

第三版

免疫原:Fusion protein of USP5

亚型:IgG2a

克隆:356

**描述:**USP5,also named as ISOT,belongs to the peptidase C19 family. Knock-down of USP5 causes the accumulation of TP53/p53 and an increase in TP53/p53 transcriptional activity because the unanchored polyubiquitin that accumulates is able to compete with ubiquitinated TP53/p53 but not with MDM2 for proteasomal recognition. USP5 is a potential target for p53 activating therapeutic agents for the treatment of cancer. It is a novel proteasome associated protein. USP5 cleaves linear and branched multiubiquitin polymers with a marked preference for branched polymers.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF