

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GC Monoclonal Antibody

单克隆抗体

产品货号: MA02284

第三版

免疫原:Fusion protein of Vitamin D binding protein

亚型:IgG1

克隆:146

描述:Vitamin D binding protein is a sparsely glycosylated serum protein responsible for highly specific binding and tissue-specific delivery of vitamin D and its metabolites. In addition, it is also an actin scavenger, and is the precursor to the immunomodulatory protein, Gc-MAF. Vitamin D binding protein has been proposed to have significant roles in C5a chemotaxis, osteoclast development and possibly in macrophage activation/recruitment.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF