



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human 14-3-3 sigma Monoclonal Antibody

单克隆抗体

产品货号: MA02285

第三版

免疫原:Fusion protein of 14-3-3 sigma

亚型:IgG1

克隆:954

描述:14-3-3 sigma is a member 14-3-3 family of highly conserved proteins participating in diverse cellular processes including cell cycle progression. It is exclusively present in various epithelial cells, and alternatively named as human epithelial marker (HEM), or stratifin. Alteration of its expression has been linked to cancer.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF