



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LIMA1 Monoclonal Antibody

单克隆抗体

产品货号: MA02290

第三版

免疫原:Fusion protein of EPLIN

亚型:IgG1

克隆:6578

描述:EPLIN is a cytoskeletal protein that is preferentially expressed in epithelial cells and is implicated in regulation of cytoskeletal dynamics and cytokinesis. In human,there are two known isoforms,EPLIN-a and -b,generated by alternative promoter usage from a single gene. EPLIN-a expression is often down-regulated in cancerous cells and tissues. This antibody recognizes both of isoforms of EPLN.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF