

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse HSD17B10 Monoclonal Antibody

单克隆抗体

产品货号: MA02291

第三版

免疫原:Fusion protein of ERAB

亚型:IgG1

克隆:9546

描述:HSD17B10 (3-hydroxyacyl-CoA dehydrogenase type-2) is a multifunctional mitochondrial enzyme that acts on a wide spectrum of substrates,including neuroactive steroids,alcohols,leucine,and fatty acids,with a preference for short-chain methyl-branched acyl-CoAs.It has 2 isoforms produced by alternative splicing.Defects in HSD17B10 are the cause of 2-methyl-3-hydroxybutyryl-CoA dehydrogenase deficiency (MHBD deficiency) and mental retardation syndromic X-linked type 10 (MRXS10).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF