



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL17D Monoclonal Antibody

单克隆抗体

产品货号: MA02295

第三版

免疫原:Fusion protein of IL-17D

亚型:IgG2b

克隆:3447

描述:IL17D is a secreted cytokine with homology to the IL-17 family of proteins. IL17D is preferentially expressed in skeletal muscle,brain,adipose tissue,heart,lung,and pancreas. The treatment of endothelial cells with IL17D has been shown to stimulate the production of other cytokines including IL6,IL8 and CSF2/ GM-CSF. The increased expression of IL8 induced by IL17D was found to be NF-kappa B-dependent.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF