
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL27 Monoclonal Antibody

单克隆抗体

产品货号: MA02297

第三版

免疫原:Fusion protein of IL-27

亚型:IgG2b

克隆:1354

描述:IL-27,a member of IL-12/IL-23 heterodimeric family of cytokines,has pleiotropic properties that can enhance or limit immune responses. IL-27 is a heterodimeric cytokine composed of two subunits: IL-27a and IL-27b. IL-27 acts on various cell types,including T cells,B cells,macrophages,dendritic cells,natural killer (NK) cells and non-hematopoietic cells. IL-27 plays a critical role in the early regulation of T helper type 1 initiation,and enhances proliferation of naive CD4+T cells and naive B cells. It,however,also exerts anti-inflammatory functions by inhibiting the development of Th17 cells and inducing IL-10 producing type 1 regulatory T cells. IL-27 is a potentially promising cytokine for therapeutic approaches on various human diseases.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF