

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ITGA3 Monoclonal Antibody

单克隆抗体

产品货号: MA02298

第三版

免疫原:Fusion protein of Integrin alpha-3

亚型:IgG2a

克隆:1269

描述:ITGA3,also named as MSK18,CD49c,GAPB3 and FRP-2,belongs to the integrin alpha chain family. It is a receptor for fibronectin,laminin,collagen,epiligrin,thrombospondin and CSPG4. ITGA3 may mediate with LGALS3 the stimulation by CSPG4 of endothelial cells migration. ITGA3 is highly concentrated in epithelial cells where it strongly adheres to Laminin-5 and Laminin-5 induced rapid adhesion can be blocked by antibodies against the alpha-3 integrin subunit. It can be cleaved into heavy chain(840aa) and light chain(176aa). This antibody can recognize the heavy chain.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF