

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLC45A3 Monoclonal Antibody

单克隆抗体

产品货号: MA02301

第三版

免疫原:Fusion protein of Prostein

亚型:IgG2b

克隆:245

描述:Prostein is a prostate tissue-specific protein whose expression is highly restricted to normal and cancerous prostate tissues. Anti-prostein is useful for the identification of prostate adenocarcinomas and is a good marker to identify prostatic origin in metastatic cancer.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF