

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat CUL4B Monoclonal Antibody

单克隆抗体

产品货号: MA02305

第三版

免疫原:Fusion protein of CUL4B

亚型:IgG1

克隆:9603

**描述:**Cullin-RING ligases (CRLs) complexes participate in the regulation of diverse cellular processes,including cell cycle progression,transcription,signal transduction and development . Serving as the scaffold protein,cullins are crucial for the assembly of ligase complexes,which recognize and target various substrates for proteosomal degradation. Two cullin 4 (CUL4) proteins,CUL4A (87 kDa) and CUL4B(104 kDa),are two members in cullin family with 83% of identity. Mutations in human CUL4B are one of the major causes of X-linked mental retardation. Cul4b knockout mice demonstrated that CUL4B is indispensable for embryonic development in the mouse .

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF