



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EPHA1 Monoclonal Antibody

单克隆抗体

产品货号: MA02307

第三版

免疫原:Synthetic Peptide

亚型:IgG1

克隆:650

描述:EPHA1,also named as EPH,EPHT and EPHT1,belongs to ephrin receptor subfamily of the protein-tyrosine kinase family. EPH and EPH-related receptors have been implicated in mediating developmental events,particularly in the nervous system. Receptors in the EPH subfamily typically have a single kinase domain and an extracellular region containing a Cys-rich domain and 2 fibronectin type III repeats. Ligands in the ephrin-A subclass,including the prototype family member ephrin-A1 (B61),are membrane associated through glycosylphosphatidyl-inositol linkages,whereas ephrin-B subclass consists of ligands with transmembrane domains. The general role of the Eph family is in mediating repulsive cell-cell interaction.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF