

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse GATA1 Monoclonal Antibody

单克隆抗体

产品货号: MA02308

第三版

免疫原:Fusion protein of GATA1

亚型:IgG1

克隆:1245

描述:GATA1 is first identified as a member of the GATA transcription factor family,whose members bind the consensus (WGATAR) binding motif . GATA1,a zinc finger DNA-binding transcription factor,plays a critical role in the normal development of hematopoietic cell lineages. The protein contains an N-terminal region that confers transcriptional activity and a C-terminal domain that mediates binding to DNA and other factors . GATA-1 is also implicated in regulating the expression of the erythroid and megakaryocytic-specific genes .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF