

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PDGFRA Monoclonal Antibody

单克隆抗体

产品货号: MA02310

第三版

免疫原:Fusion protein of PDGFRA

亚型:IgG1

克隆:134

描述:Platelet-derived growth factor receptor alpha (PDGFRA,also known as CD140a) is a cell surface tyrosine kinase receptor for members of the platelet-derived growth factor family. PDGFRA is implicated in the regulation of embryonic development,tumorigenesis,cell proliferation,survival and chemotaxis. Mutations in PDGFRA gene have been associated with idiopathic hypereosinophilic syndrome,somatic and familial gastrointestinal stromal tumors,and a variety of other cancers.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF