



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse UBE3A Monoclonal Antibody

单克隆抗体

产品货号: MA02312

第三版

免疫原:Fusion protein of UBE3A

亚型:IgG1

克隆:756

描述:UBE3A(Ubiquitin-protein ligase E3A) is also named as E6AP,EPVE6AP,HPVE6A and belongs to the E3 ubiquitin protein ligase family. It functions as both an E3 ligase in the ubiquitin proteasome pathway and as a transcriptional coactivator. It is also initially identified as a cellular protein that mediates in vitro association of the human papillomavirus E6 protein with p53,leading to the ubiquitin-dependent degradation of p53. Defects in UBE3A are a cause of Angelman syndrome (AS). It has 3 isoforms produced by alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF