

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

VSV-G-Tag Monoclonal Antibody

单克隆抗体

产品货号: MA02320

第三版

免疫原:Synthetic Peptide

亚型:IgG

克隆:1H7

描述: Vesicular stomatitis virus (VSV), an enveloped RNA virus from the Rhabdoviridae family, is released from the plasma membrane of host cells by a process called budding. The glycoprotein (VSV-g) contains a domain in its extracellular membrane proximal stem that appears to be needed for efficient VSV budding. This antibody can be used to detect the VSVG tagged protein.

配方:

浓度:1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF