



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

GST-Tag Monoclonal Antibody

单克隆抗体

产品货号: MA02321

第三版

免疫原:Recombinant Protein

亚型:IgG

克隆:2B7

**描述:**GST(Glutathione S-Transferase) is a widely used protein tag encoded by the Schistosoma japonicum. GST provides both an easily detectable tag and a simple purification process with little effect on the biological function of the protein of interest. This GST-fusion protein can then be purified from cells via its high affinity for glutathione. Antibodies to GST are useful for checking protein expression both in plaques and on western blots as well as for immunoaffinity purification of proteins expressed in cells. This antibody can be used to detect the GST tagged protein.

**配方:**

**浓度:**1mg/mL

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF