



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LE/AF Purified Anti-Human CD4 Monoclonal Antibody

单克隆抗体

产品货号: MA02375

第三版

免疫原:None

亚型:Mouse IgG2b

克隆:OKT4

描述:CD4, also known as T4, is a 55 kD single-chain type I transmembrane glycoprotein expressed on most thymocytes, a subset of T cells, and monocytes/macrophages. CD4, a member of the Ig superfamily, recognizes antigens associated with MHC class II molecules and participates in cell-cell interactions, thymic differentiation, and signal transduction. CD4 acts as a primary receptor for HIV, binding to HIV gp120. CD4 has also been shown to interact with IL-16. The OKT4 antibody binds to the D3 domain of CD4 and does not block HIV binding.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF