



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LE/AF Purified Anti-Human CD8 Monoclonal Antibody

单克隆抗体

产品货号: MA02376

第三版

免疫原:None

亚型:Mouse IgG2a

克隆:OKT8

描述:CD8 is known as T8 and Leu2, is a member of the immunoglobulin superfamily found on the majority of thymocytes, a subset of peripheral blood T cells, and NK cells (which express almost exclusively CD8a homodimers). CD8 acts as a co-receptor with MHC class I-restricted T cell receptors in antigen recognition and T cell activation and has been shown to play a role in thymic differentiation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF