



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LE/AF Purified Anti-Human CD40 Monoclonal Antibody

单克隆抗体

产品货号: MA02378

第三版

免疫原:None

亚型:Mouse IgG1

克隆:G28.5

**描述:**CD40 is a 48 kD type I glycoprotein also known as BP50. It is a member of the TNFR superfamily primarily expressed on B cells, macrophages, follicular dendritic cells, endothelial cells, fibroblasts, and at low levels on plasma cells. CD40 has been reported to be involved in B cell differentiation, costimulation, isotype class-switching, and protection from apoptosis. Additionally, CD40 is important for T cell-B cell interactions. The ligand of CD40 is CD154 (CD40 ligand). The G28.5 antibody has been reported to promote B cell proliferation in the presence of anti-IgM, IL-4 or PMA.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF