



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Purified Anti-Human CD10 Monoclonal Antibody

单克隆抗体

产品货号: MA02383

第三版

免疫原:None

亚型:Mouse IgG1

克隆:HI10a

**描述:**CD10 is a 100 kD neutral endopeptidase and a member of the metalloprotease family. It is a type II transmembrane protein also known as common acute lymphoblastic leukemia antigen (CALLA), enkephalinase, and neprilysin. CD10 is expressed on B cell precursors, T cell precursors, and neutrophils. CD10 is involved in B cell development and has been shown to bind opioid enkephalins, bradykinin, angiotensins I and II, and other biologically active peptides.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF