



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Purified Anti-Human CD19 Monoclonal Antibody

单克隆抗体

产品货号: MA02385

第三版

免疫原:None

亚型:Mouse IgG1

克隆:HIB19a

**描述:**CD19 is a 95 kD type I transmembrane glycoprotein is also known as B4. It is a member of the immunoglobulin superfamily expressed on B-cells (from pro-B to blastoid B cells, absent on plasma cells) and follicular dendritic cells. CD19 is involved in B cell development, activation, and differentiation. CD19 forms a complex with CD21 (CR2) and CD81 (TAPA-1), and functions as a BCR co-receptor.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF