



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Purified Anti-Human CD34 Monoclonal Antibody

单克隆抗体

产品货号: MA02387

第三版

免疫原:None

亚型:Mouse IgG2b

克隆:4H11

**描述:**CD34, also known as gp105-120, is a type I monomeric sialomucin-like glycoprophoprotein with an approximate molecular weight of 105-120 kD. It is selectively expressed on the majority of hematopoietic stem/progenitor cells, bone marrow stromal cells, capillary endothelial cells, embryonic fibroblasts, and some nervous tissue. CD34 is a commonly used marker for identifying human hematopoietic stem/progenitor cells. Based on different sensitivities, four groups of epitopes of CD34 have been described. CD34 mediates cell adhesion and lymphocytes homing through binding to L-selectin and E-selectin ligands.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF