



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Purified Anti-Human CD56 Monoclonal Antibody

单克隆抗体

产品货号: MA02392

第三版

免疫原:None

亚型:Mouse IgG1

克隆:B-A19

描述:CD56 is a single transmembrane glycoprotein also known as NCAM (Neural Cell Adhesion Molecule), Leu-19, or NKH1. The CD56 (B-A19) antibody recognizes a 175-220 kDa protein, the neural cell adhesion molecule 1 (NCAM-1), a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation and neurite outgrowth.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF