

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse LE/AF Purified Anti-Mouse CD16/CD32 Monoclonal Antibody

单克隆抗体

产品货号: MA02399

第三版

免疫原:None

亚型:Rat IgG2b

克隆:2.4G2

描述:CD16 is the low affinity IgG Fc receptor III (FcR III) and CD32 is FcR II. CD16/CD32 are expressed on B cells, monocytes/macrophages, NK cells, granulocytes, mast cells, and dendritic cells. The Fc receptors bind antibody-antigen immune complexes and mediate adaptive immune responses. This 2.4G2 mAb is useful for blocking non-specific binding of immunoglobulin to the Fc receptors.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF