



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse LE/AF Purified Anti-Mouse CD90 Monoclonal Antibody

单克隆抗体

产品货号: MA02407

第三版

免疫原:None

亚型:Rat IgG2a

克隆:M5/49.4.1

描述:CD90 is a 25-35 kD immunoglobulin superfamily member, also known as Thy-1. It is expressed on hematopoietic stem cells and neurons, all thymocytes, and peripheral T cells, but not on B cells. CD90 is a glycosylphosphatidylinositol (GPI)-anchored membrane glycoprotein involved in signal transduction. CD90 is involved in costimulation of lymphocyte proliferation and activation, and hematopoietic stem cell differentiation. CD90 has been shown to interact with CD45.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF