



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse LE/AF Purified Anti-Mouse CD137 Monoclonal Antibody

单克隆抗体

产品货号: MA02409

第三版

免疫原:None

亚型:Rat IgG2a

克隆:LOB12.3

**描述:**CDw137 is a 39 kD TNF-receptor superfamily member also known as 4-1BB, TNFRSF9, Ly-63, and ILA (Induced by Lymphocyte Activation). CDw137 is expressed as monomers, dimers, and tetramers on activated T and B cells. Interaction of CDw137 with its ligand, 4-1BBL (CDw137L), is important in T-B cell costimulation. The 17B5 antibody has been reported to block 4-1BBL costimulated T cell proliferation.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF