



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse LE/AF Purified Anti-Mouse CD137L Monoclonal Antibody

单克隆抗体

产品货号: MA02410

第三版

免疫原:None

亚型:Rat IgG2a

克隆:TKS-1

描述:4-1BB ligand, also known as CD137L, is a 97 kD member of the TNF superfamily present on antigen-presenting cells (dendritic cells and macrophages) and activated B and T lymphocytes. 4-1BBL interacts with 4-1BB (CDw137) to induce T and B cell costimulation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF