



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse LE/AF Purified Anti-Mouse CD154 Monoclonal Antibody

单克隆抗体

产品货号: MA02412

第三版

免疫原:None

亚型:Hamster IgG

克隆:MR1

描述:CD154 is a 39 kD TNF superfamily member also known as CD40 ligand, gp39, T-BAM, TRAP, and Ly-62. CD154 is an accessory molecule expressed predominantly on activated CD4+ lymphocytes that bind CD40. CD154 plays an important role in T-B cell costimulation. The MR1 antibody has been reported to inhibit the activation of T and B lymphocytes in vitro and antigen-specific lymphocyte responses in vivo.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF