



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse LE/AF Purified Anti-Mouse IL-4 Monoclonal Antibody

单克隆抗体

产品货号: MA02419

第三版

免疫原:None

亚型:Rat IgG1

克隆:11B11

描述:IL-4 is a pleiotropic cytokine produced by activated T cells, mast cells, and basophils. IL-4 is a potent lymphoid cell growth factor which stimulates the growth and activation of certain B cells and T cells. IL-4 is important for regulation of T helper subset development.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF