

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse/Human LE/AF Purified Anti-Mouse/Human IL-7 Monoclonal Antibody

单克隆抗体

产品货号: MA02421

第三版

免疫原:None

亚型:Mouse IgG2b

克隆:M25

**描述:**The M25 monoclonal antibody reacts with mouse and human IL-7, a 20-28 kDa hematopoietic growth factor secreted by stromal cells in the bone marrow and thymus. IL-7 is also produced by keratinocytes, dendritic cells, hepatocytes, neurons, and epithelial cells. IL-7 binds to the IL-7 receptor, a heterodimer consisting of IL-7 receptor alpha and the common gamma chain receptor. Binding results in stimulation of differentiation of pluripotent hematopoietic stem cells into lymphoid progenitor cells and proliferation of all lymphoid cells. IL-7 is also important for proliferation of pre-B cells as well as T and NK cell survival, development and homeostasis. The M25 antibody is useful for neutralizing the bioactivity of IL-7 as well as administration of IL-7 cytokine. IL-7 administration as a complex with M25 prolongs IL-7 availability in vivo by decreasing specific and nonspecific consumption.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF